

QuickVue At-Home COVID-19 Test B-Roll Package Video – Descriptive Transcript

Visual (no audio)
Quidel logo animates in.
Transition to the 2-Count QuickVue At-Home COVID-19 Test product package, a prefilled tube, nasal swab, test strip and timer on a counter.
Transition to a woman's hand removing the prefilled tube, nasal swab, and concealed test strip from the tray.
The woman picks up a prefilled tube and twists the cap to open the tube. The woman placed the cap on the counter and then places the tube in a hole within the tray.
The woman peels open the swab wrapper and pulls out the swab by its handle.
The woman inserts the swab and swirls it inside her right nostril.
The woman inserts the swab into the prefilled tube and swirls.
The woman inserts the test strip into the tube, with the arrows pointing downward into the tube.
The woman places the test strip under Step 10 within the outline included on the User Instructions sheet.
Transition to the package, test kit components and timer on the counter.
Quidel logo animates in.
Copyright date appears at the bottom of the screen “© 2008-2021 Quidel Corporation”